

L-753, EHRlich'S HEMATOXYLIN

Fixation: Formalin or fixative of choice.

Sections: Paraffin

Staining:

1. Deparaffinize and hydrate to distilled water.
2. Stain for 2-5 minutes in **Ehrlich's Hematoxylin** (L-753-1). Wash in water.
3. Differentiate in **Acid Alcohol** (L-753-2) until the sections turn red.
4. Blue sections in **Sodium Acetate, 1%** (L-753-3) or **Ammonia Water** (L-753-3A) or **Lithium Carbonate Solution** (L-753-3B).
5. Wash thoroughly and counterstain.

REFERENCES:

Ehrlich, P., Biol. Centralbl., 6:214:224, 1886

Mallory, F.B., Pathological Technique, W.B. Saunders Co. Philadelphia, p 73, 1938

Clark, G., (ed.) Staining Procedures, 3rd edition, Williams & Wilkins, Blatimore, pp 33-34, c 1973

Luna, L.G., (Ed.), Manual of Histologic Staining Methods of the AFIP, 3rd Edition, McGraw-Hill, NY, pp 35 and 38, c 1968.