




p53

USA IVD

This recognizes a 53kDa phosphoprotein. Anti P53 positivity has also been used to differentiate uterine serous carcinoma from endometrioid carcinoma.

Mouse Monoclonal, Clone: D07

Control: Colon Carcinoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG2b/K **Cellular localization:** Nuclear



Colon Carcinoma

Catalog no.

PF-2010-1

PF-2010-2

PF-2010-3

PF-2010-4

PF-2010-5

PF-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides