




# A-1-Antichymotrypsin

USA IVD

Its major application is defining the presence of Alpha-1-Antichymotrypsin in histiocytes and tumors derived from them

## Rabbit Polyclonal

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Isotype:** **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.	Volume
A1A-2010-1	0.1 ml concentrated
A1A-2010-2	0.5 ml concentrated
A1A-2010-3	1 ml concentrated
A1A-2010-4	1 ml Prediluted
A1A-2010-5	7 ml Prediluted
A1A-2010	5 Positive control slides