




Androgen Receptor

USA IVD

This has been useful in differentiating morpheaform basal cell carcinoma from desmoplastic trichoepithelioma in the skin.

Mouse Monoclonal, Clone: AR441

Control: Prostate. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG1 **Cellular localization:** Nuclear



Prostate

| Catalog no. | Volume |
|-------------|---------------------------|
| AR-2010-1 | 0.1 ml concentrated |
| AR-2010-2 | 0.5 ml concentrated |
| AR-2010-3 | 1 ml concentrated |
| AR-2010-4 | 1 ml Prediluted |
| AR-2010-5 | 7 ml Prediluted |
| AR-2010 | 5 Positive control slides |