




## CD105 (Endoglin)

USA IVD

This could serve as a target for anti-angiogenesis therapy.

**Mouse Monoclonal. Clone: MRQ-14**

**Control:** Renal Cell Carcinoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.  
**Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K. **Cellular localization:** Cytoplasmic



Renal Cell Carcinoma

Catalog no.	Volume
C5E-2010-1	0.1 ml concentrated
C5E-2010-2	0.5 ml concentrated
C5E-2010-3	1ml concentrated
C5E-2010-4	1ml Prediluted
C5E-2010-5	7 ml Prediluted
C5E-2010	5 Positive control slides