




## CD2

USA IVD

This is useful for the identification of all normal T-lymphocytes.

**Mouse Monoclonal, Clone: MRQ-11**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic , Membranous.



Tonsil

Catalog no.	Volume
CD2-2010-1	
CD2-2010-2	0.1 ml concentrated
CD2-2010-3	0.5 ml concentrated
CD2-2010-4	1 ml concentrated
CD2-2010-5	1 ml Prediluted
CD2-2010	7 ml Prediluted
	5 Positive control slides

