




## CD4

USA IVD

This is used in the immunophenotyping of reactive lymphocytes and lymphoproliferative disorders.

**Rabbit Monoclonal, Clone: SP35**

**Control:** Tonsil . **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Membranous, Cytoplasmic



Tonsil

<b>Catalog no.</b>	<b>Volume</b>
CD4-2010-1	0.1 ml concentrated
CD4-2010-2	0.5 ml concentrated
CD4-2010-3	1 ml concentrated
CD4-2010-4	1 ml Prediluted
CD4-2010-5	7 ml Prediluted
CD4-2010	5 Positive control slides

