



CD45RO

USA IVD

This is a good marker for T-cell tumors because it shows no reactivity with B-cells.

Mouse Monoclonal, Clone: UCHL-1

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2a}/K Cellular **localization:** Membranous



Tonsil

Catalog no.	Volume
45O-2010-1	0.1 ml concentrated
45O-2010-2	0.5 ml concentrated
45O-2010-3	1ml concentrated
45O-2010-4	1ml Prediluted
45O-2010-5	7 ml Prediluted
45O-2010	5 Positive control slides