




CD99

USA IVD

This can be found on the cell membrane of Ewing's Sarcoma and primitive peripheral neuroectodermal tumors.

Mouse Monoclonal. Clone: H036-1.1

Control: Ewing's Sarcoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** H036-1.1 IgM. **Cellular localization:** Cytoplasmic, Membranous



Ewing's Sarcoma

Catalog no.	Volume
C99-2010-1	0.1 ml concentrated
C99-2010-2	0.5 ml concentrated
C99-2010-3	1ml concentrated
C99-2010-4	1ml Prediluted
C99-2010-5	7 ml Prediluted
C99-2010	5 Positive control slides