




# Collagen Type IV

USA IVD

This can be useful in the classification of soft tissue tumors: schwannomas and leiomyomas.

**Mouse Monoclonal, Clone: CIV22**

**Control:** Lung. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Intercellular



Lung

Catalog no.	Volume
Col-2010-1	0.1 ml concentrated
Col-2010-2	0.5 ml concentrated
Col-2010-3	1ml concentrated
Col-2010-4	1ml Prediluted
Col-2010-5	7 ml Prediluted
Col-2010	5 Positive control slides