




Desmin

USA IVD

This is used to demonstrate the myogenic components and derivation of tumors.

Mouse Monoclonal, Clone: D33

Control: Muscle. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Muscle

| Catalog no. | Volume |
|-------------|---------------------------|
| De-2010-1 | 0.1 ml concentrated |
| De-2010-2 | 0.5 ml concentrated |
| De-2010-3 | 1ml concentrated |
| De-2010-4 | 1ml Prediluted |
| De-2010-5 | 7 ml Prediluted |
| De-2010 | 5 Positive control slides |

