



Epithelial Specific Antigen/Ep-CAM (MOC-

31)

USA IVD

This is used to distinguish adenocarcinoma from mesothelioma and hepatocellular carcinoma.

Mouse Monoclonal, Clone: MOC-31

Control: Colon Carcinoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Colon Carcinoma

Catalog no.

ES-2010-1

Volume

ES-2010-2

0.1 ml concentrated

ES-2010-3

0.5 ml concentrated

ES-2010-4

1ml concentrated

ES-2010-5

1ml Prediluted

ES-2010

7 ml Prediluted

5 Positive control slides

