




# FOXP1

USA IVD

This may be a diagnostic marker of B-cell lymphoma, leg type.

## Rabbit Polyclonal

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Cellular localization:** Nuclear



Tonsil

<b>Catalog no.</b>	<b>Volume</b>
FX-2010-1	0.1 ml concentrated
FX-2010-2	0.5 ml concentrated
FX-2010-3	1ml concentrated
FX-2010-4	1ml Prediluted
FX-2010-5	7 ml Prediluted
FX-2010	5 Positive control slides