




## Factor XIIIa

USA IVD

This is used to differentiate between dermatofibroma, dermafibrosarcoma protuberans and desmoplastic malignant melanoma.

**Rabbit Monoclonal, Clone: EP3372**

**Control:** Dermatofibroma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** AC-1A1, IgG<sub>1</sub>, EP3372, IgG **Cellular localization:** Cytoplasmic



Dermatofibroma

<b>Catalog no.</b>	<b>Volume</b>
Fr-2010-1	0.1 ml concentrated
Fr-2010-2	0.5 ml concentrated
Fr-2010-3	1ml concentrated
Fr-2010-4	1ml Prediluted
Fr-2010-5	7 ml Prediluted
Fr-2010	5 Positive control slides

