




Follicular Dendritic Cell

USA IVD

This is useful in the identification of follicular dendritic cell matrix found in normal lymph nodes and tonsil tissue.

Mouse Monoclonal, Clone: CNA.42

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgM/K **Cellular localization:** Cytoplasmic



Tonsil

| Catalog no. | Volume |
|--------------------|---------------------------|
| FD-2010-1 | 0.1 ml concentrated |
| FD-2010-2 | 0.5 ml concentrated |
| FD-2010-3 | 1ml concentrated |
| FD-2010-4 | 1ml Prediluted |
| FD-2010-5 | 7 ml Prediluted |
| FD-2010 | 5 Positive control slides |