




## Galectin-3

USA IVD

This used for the differentiation between follicular adenoma, follicular carcinoma, and papillary carcinoma of thyroid.

**Mouse Monoclonal: 9C4**

**Control:** Papillary Carcinoma of Thyroid. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic, Nuclear



Papillary Carcinoma of Thyroid

Catalog no.	Volume
GA-2010-1	0.1 ml concentrated
GA-2010-2	0.5 ml concentrated
GA-2010-3	1ml concentrated
GA-2010-4	1ml Prediluted
GA-2010-5	7 ml Prediluted
GA-2010	5 Positive control slides