




## Granzyme B

USA IVD

This is useful in diagnosing natural killer/T-cell lymphoma, and anaplastic large cell lymphoma.

### Rabbit Polyclonal

**Control:** Lymphoblastic Lymphoma, Spleen **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Granular Cytoplasmic



Lymphoblastic Lymphoma, Spleen

### Catalog no.

GB-2010-1

GB-2010-2

GB-2010-3

GB-2010-4

GB-2010-5

GB-2010

### Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides