




# HHV-8

USA IVD

This labels the latent nuclear antigen protein via immunohistochemistry. It have been found in Kaposi's sarcoma.

**Mouse Monoclonal, Clone: 13B10**

**Control:** Kaposi's Sarcoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Nuclear



Kaposi's Sarcoma

## Catalog no.

HH-2010-1

HH-2010-2

HH-2010-3

HH-2010-4

HH-2010-5

HH-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .