



Hepatitis B Virus Surface Antigen


USA ASR

This labels the cell cytoplasm infected with hepatitis B virus.

Mouse Monoclonal, Clone: S1-210

Control: Hepatitis B Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

Shelf life: Up to 36 months. **Isotype:** IgG_{2b} **Cellular localization:** Cytoplasmic.



Hepatitis B Infected Tissue

Catalog no.

Volume

HE-2010-1	0.1 ml concentrated
HE-2010-2	0.5 ml concentrated
HE-2010-3	1ml concentrated
HE-2010-4	1ml Prediluted
HE-2010-5	7 ml Prediluted
HE-2010-6	5 Positive control slides.