




Hepatocyte Specific Antigen

USA IVD

This is useful in differentiating hepatocellular carcinomas with adenoid features from adenocarcinomas.

Mouse Monoclonal, Clone: OCH1E5

Control: Liver **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Liver

Catalog no.

HS-2010-1

Volume

HS-2010-2

0.1 ml concentrated

HS-2010-3

0.5 ml concentrated

HS-2010-4

1ml concentrated

HS-2010-5

1ml Prediluted

HS-2010

7 ml Prediluted

5 Positive control slides.