




Insulin

USA IVD

When staining the cytoplasm of tumors, this is the most reliable indication of functional insulinomas.

Guinea Pig Polyclonal

Control: Pancreas **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Pancreas

Catalog no.

IS-2010-1

IS-2010-2

IS-2010-3

IS-2010-4

IS-2010-5

IS-2010

1ml Prediluted

7 ml Prediluted

5 Positive control slides .

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated