




MUM1

USA IVD

This is a good marker of Hodgkin's and Reed-Sternberg cells of classic Hodgkin's disease.

Mouse Monoclonal, Clone: MRQ-8

Control: Tonsil Embedded: Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic, Nuclear



Tonsil

Catalog no.

MM-2010-1

MM-2010-2

MM-2010-3

MM-2010-4

MM-2010-5

MM-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.