




Macrophage (HAM-56)

USA IVD

This interacts with tangible macrophages, interdigitating macrophages of lymph nodes, and tissue macrophages.

Mouse Monoclonal: HAM-56

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgM/K **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.

MA-2010-1

MA-2010-2

MA-2010-3

MA-2010-4

MA-2010-5

MA-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .