




# Myeloperoxidase

USA IVD

This detects granulocytes in blood and precursors of granulocytes in bone marrow.

## Rabbit Polyclonal

**Control:** Bone Marrow **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Bone Marrow

## Catalog no.

MY-2010-1

MY-2010-2

MY-2010-3

MY-2010-4

MY-2010-5

MY-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .