




PD-1

USA IVD

This is a marker for angioimmunoblastic lymphoma.

Mouse Monoclonal, Clone: MRQ-22

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.

PD-2010-1

PD-2010-2

PD-2010-3

PD-2010-4

PD-2010-5

PD-2010-6

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides