




PGP 9.5

USA IVD

Rabbit Polyclonal

Control: Nerve Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Nerve Tissue

Catalog no.

PG-2010-1

PG-2010-2

PG-2010-3

PG-2010-4

PG-2010-5

PG-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides .