



# PSA

USA IV D

This recognizes primary and metastatic prostatic neoplasms and rarely tumors of the neoplastic origin.

**Mouse Monoclonal, Clone: ER-PR8**

**Control:** Prostate **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic



Prostate

## Catalog no.

PA-2010-1

PA-2010-2

PA-2010-3

PA-2010-4

PA-2010-5

PA-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides