




PU.1

USA IVD

This is a transcription factor that is important for B-cell development. Various lymphomas are also positive for this marker.

Rabbit Monoclonal, Clone: EPR3158Y

Control: Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Nuclear



Tonsil

Catalog no.

PU-2010-1

PU-2010-2

PU-2010-3

PU-2010-4

PU-2010-5

PU-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides