




# Perforin

USA IVD

This is a pore-forming protein that leads to osmotic lysis of the target cells and enables granzymes to enter the target cells and activate apoptosis.

**Mouse Monoclonal, Clone: MRQ-23**

**Control:** Spleen **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Perinuclear, Cytoplasmic



Spleen

## Catalog no.

PE-2010-1

PE-2010-2

PE-2010-3

PE-2010-4

PE-2010-5

PE-2010

## Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides