




## Podoplanin (D2-40)

USA IVD

This is a transmembrane mucoprotein. Podoplanin is selectively expressed in lymphatic endothelium as well as lymphangiomas Kaposi sarcomas and in a subset of angiosarcomas.

**Mouse Monoclonal, Clone: C2-40**

**Control:** Tonsil **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Tonsil

### Catalog no.

PN-2010-1

PN-2010-2

PN-2010-3

PN-2010-4

PN-2010-5

PN-2010

### Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides