




Progesterone Receptor

USA ASR

This is used as a prognostic indicator in breast carcinoma.

Rabbit Monoclonal, Clone: Y85

Control: Breast **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Nuclear



Breast

Catalog no.

PR-2010-1

PR-2010-2

PR-2010-3

PR-2010-4

PR-2010-5

PR-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides