



SV40


USA ASR

This is associated with brain tumors, bone cancers, non-Hodgkin's lymphoma, and malignant mesothelioma.

Mouse Monoclonal, Clone: MRQ-4

Control: SV40 Infected Tissue **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

Shelf life: Up to 36 months. **Isotype:** IgG_{2a} **Cellular localization:** Nuclear



SV40 Infected Tissue

Catalog no.

SV-2010-1

SV-2010-2

SV-2010-3

SV-2010-4

SV-2010-5

SV-2010

Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides