



# TCL 1


USA IV D

This is useful with the discrimination of B-cell lymphomas from T-cell lymphomas, CD30+ anaplastic large cell lymphomas, multiple myeloma, and marginal zone B-cell lymphoma.

**Mouse Monoclonal, Clone: MRQ-7**

**Control:** Tonsil, B-cell Lymphoma **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Cytoplasmic, Nuclear



Tonsil

## Catalog no.

TC-2010-1

TC-2010-2

TC-2010-3

TC-2010-4

TC-2010-5

TC-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides.