




# Tryptase

USA IV D

This is a marker for mast cells, basophils, and their derivatives.

**Mouse Monoclonal, Clone: G3**

**Control:** Uterus **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic



Uterus

## Catalog no.

TP-210-1

TP-2010-2

TP-2010-3

TP-2010-4

TP-2010-5

TP-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides