




# Uroplakin III

USA IV D

This is a better indicator of the malignant potential of the tumor than the grade of the tumor.

**Mouse Monoclonal, Clone: AU-1**

**Control:** Bladder **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic, Membranous



Bladder

## Catalog no.

UR-210-1

UR-2010-2

UR-2010-3

UR-2010-4

UR-2010-5

UR-2010

## Volume

0.1 ml Concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides