




## bcl-6

USA IVD

This is a transcriptional regulator gene which codes for a 706-amino-acid nuclear zinc finger protein.

**Mouse Monoclonal, Clone: GI191E/A8**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Nuclear



Tonsil

<b>Catalog no.</b>	<b>Volume</b>
BC6-210-1	0.1 ml concentrated
BC6-210-2	0.5 ml concentrated
BC6-210-3	1 ml concentrated
BC6-210-4	1 ml Prediluted
BC6-210-5	7 ml Prediluted
BC6-210	5 Positive control slides

