

DETECTION SYSTEMS

Product	Size	Cat.#
RBI Master Detection™ HRP Polymer System	50ml (500 tests)	RB-68-1
	100ml (1,000 tests)	RB-68-2
RBI Master Detection™ Alk-Phos.Polymer system	50 ml (500 tests)	RB-69-1
	100 ml (1,000 tests)	RB-69-2
RBI-Completed™ HRP Polymer Complete Detection system with DAB chromagen.	15 ml (150 tests)	RB-70-1
	50 ml (500 test)	RB-70-2
	100ml (1,000 tests)	RB-70-3
HRP Avidin Biotin	50 ml (500 tests)	RB-71-1
	200 ml (2,000 tests)	RB-71-2
Alk. Phos. Avidin Biotin	50 ml (500 tests)	RB-72-1
	200 ml (2,000 tests)	RB-72-2

Chromagens for HRP Systems

DAB KIT:.....	1 ml DAB/20 ml buffer(200 test).....	RB-73-1
	3 ml DAB/50 ml buffer (500 tests).....	RB-73-2
	12 ml DAB/200 ml buffer (2,000 tests).....	RB-73-3
	30 ml DAB/500 ml buffer (5,000 tests).....	RB-73-2
AEC ready-to-use:.....	50 ml (500 tests)	RB-74-1
	200 ml (2,000 tests).....	RB-74-2

Chromagens for Alk. Phos. Systems

Permanent Magenta Kit:.....	2 reagents at 1 ml ea./30 ml substrate (300 tests).....	RB-75-1
	2 reagents at 4 ml ea./100 ml substrate (1,000 tests).....	RB-75-2
Permanent Red Kit:.....	3 reagents at 0.7 ml ea./30 ml buffer (300 tests).....	RB-76-1
	3 reagents at 2.25 ml ea./100 ml buffer (300 tests).....	RB-76-2
BCIP/NBT Kit:.....	3 reagents at 0.7 ml ea./30 ml buffer (300 tests).....	RB-77-1

IHC-REAGENTS.

RBI –Master pre-treatment. DRU-2010™ (Deparaffinize, Rehydrate and Unmask.

DRU-2010.....	8 ounces.....	RBM-1
	1Pint.....	RBM-2
	1 quart.....	RBM-3
	1 gallon.....	RBM-3
IHC wash buffer.....	8 ounces.....	IHC-1
	1 pint.....	IHC-2
	1 quart.....	IHC-3
Antibody Diluent.....	50 ml.....	AD-1
	250 ml.....	AD-2
Mouse Negative Control Serum.....	15 ml.....	MN-1
Rabbit Negative Control Serum.....	15 ml.....	RN-1
Peroxide Block.....	15 ml.....	PB-1
	200ml.....	PB-2
	1 quart.....	PB-3

Background Block.....	15 ml.....	BB-1
	200ml.....	BB-2
	1 quart.....	BB-3
Avidin/Biotin Block.....	1 ml.....	AV-1
	15 ml.....	AV-2
Protease Digestive Enzyme.....	15 ml.....	PD-1
	50 ml.....	PD-2