




p63

USA RUO

This is used to differentiate benign versus malignant prostate and breast lesions.

Mouse Monoclonal, Clone: 4A4

Control: Prostate **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2a}/K **Cellular localization:** Nuclear



Prostate

Catalog no.

PS-2010-1

PS-2010-2

PS-2010-3

PS-2010-4

PS-2010-5

PS-2010

Volume

0.1 ml concentrated

0.5 ml concentrated

1ml concentrated

1ml Prediluted

7 ml Prediluted

5 Positive control slides